

# Reusable Outer Canister

The reusable out canister providing safe, convenient and economical solutions for fluid collection.

Reusable clear plastic canister used to hold the disposable semirigid liners.

Canister for wall mount 3000ml are available.

Reusable polycarbonate canister/bracket assembly for use with disposable liner.

Bracket assembly has preattached tubing for connection to disposable liner.



Product code	Description
GC1962A30	3000ml

# ABC Mouthpiece and Filter Kit

The mouthpiece and filter kit is intended to be used with Elekta Ltd. active breathing coordinator (ABC) for external beam radiation therapy treatment, each mouthpiece comes individually packaged and is a one-time-use item that should be replaced in between treatments. It comes with a filter, nose clips, mouthpiece and tubing.

Mouthpiece and connector: 22mm M

Elbow: 22mm M to 22mm F

Flex tube: 22mm F 16"

Breathing filter: 22mm M



Product code

GC8049A01

**GCMEDICA**<sup>®</sup>  
SINCE 1998

**GCMEDICA ENTERPRISE LTD. (WUXI)**

Web: <http://www.gcmmedica.com> E-mail: [info@gcmmedica.com](mailto:info@gcmmedica.com)  
Tel: 86-510-82733513 Fax: 86-51082730717  
Address: Loujin Industrial Park, Shuofang, Wuxi, 214143, China

CE 0197 FDA

# HSG Catheter

Designed for hysterosalpingography and sonohysterography procedure.  
Disposable, with rigid sheath, available with fr 5 and fr 7 in size.  
Radiopaque material of catheter allows the clear view during the procedure.



Product code	Description
GC7005A11	FR 5, latex free.
GC7005A12	FR 7, latex free.

# Suction Irrigator

Atraumatic and Non-conductive , 5mm or 10mm diameters, 350mm or 450mm lengths.  
Ergonomic handle design allows precise and single-hand control, for right or left-handed users both.  
Color coded soft PVC tubing, phthalate free.



Product code	Description
GC080741	5.4*430mm

# Rigid Suction Canister

Designed for preventing backflow and premature shut-off.

Include a patent port with 20-degree angle for 1200cc or straight port for 800cc to prevent suction tubing from kinking.

Usability assurance for collecting and retaining aspirated fluids.

Disposable, available with suction connecting tubing version.



Product code	Description
GC1957808	800cc canister, with lid
GC1957812	1200cc canister, with lid

**GC MEDICA**  
SINCE 1998

GC MEDICA ENTERPRISE LTD. IWUXII

Web: <http://www.gcmedica.com>

E-mail: [info@gcmedica.com](mailto:info@gcmedica.com)

Tel: 86-310-82733313

Fax: 86-31082733717

Address: Gujin Industrial Park, Shoufeng, Wuxi, 214163, China

CE 0187 FDA

# ENT Single Use Suction Tube

The professional range of suction offers an extensive choice of length, lumen size and tip finish to enable you to find the right match for procedure patient and surgeon.

A critical handle design provides optimum ergonomics with slim, smooth contours for comfort and improved visibility.

Engineered to prevent blocking.

Angled stainless steel cannula with suction control.

The suction control hole allows fine manipulation of suction with minimal noise.

Angled green plastic luer connector.

Automatic tip with or without side control.

Wide range of lengths and diameters are available.



Standard Suction Tip  
with Side Control



Standard Suction Tip



Olive Suction Tip



Product code	Description
OC1259A40W	FR12, Olive suction tip, tube length: 100mm
OC1259B20W	FR04, standard suction tip, tube length: 120mm
OC1259B40W	FR12, standard suction tip, tube length: 120mm
OC1259C20W	FR06, standard suction tip, tube length: 140mm
OC1259C40W	FR12, standard suction tip, tube length: 140mm
OC1259D20W	FR06, standard suction tip, tube length: 160mm
OC1259D40W	FR12, standard suction tip, tube length: 160mm
OC1259E20W	FR06, standard suction tip, tube length: 250mm, with side control at tip